

# FDA Submission

**Your Name:** Dieter Verbeemen

**Name of your Device:** The Chest X-Rays Pneumonia Notifier

## Algorithm Description

### 1. General Information

**Intended Use Statement:** Notify (help) radiologist detect pneumonia

**Indications for Use:**

- The image must be a X-ray
- The X-ray image must be a DICOM file
- The X-ray image must be taken from the chest
- The X-ray image must be taken in the AP or PA position
- The patient, man or women should be between 1 and 90 years old


**Device Limitations:**

- Any modern computer with a standard CPU, although a GPU is advised

**Clinical Impact of Performance:**

- The best use for this model is to prioritize the worklist.
- The model has a low precision and a higher recall.
  - (Recall) This means if we have 100 patients with pneumonia, Then we can properly classify these 100 patients as a positive event. (without classifying them too much as *no pneumonia*).
  - Unfortunately, the precision is not very high. This means that we will also classify healthy patients as pneumonia.
  - However, if we want to prioritize our worklist then we should prioritize our ill patients above the healthy patients (recall). But it isn't harmful to classify a healthy person with pneumonia.

### 2. Algorithm Design and Function



### DICOM Checking Steps:

- Modality must be "DX"
- Body part examined must be "CHEST"
- Patient Position must be "PA" or "AP"

### Preprocessing Steps:

- Image is normalized between 0 and 1
- Image is reshaped to a 224 x 224 image
- Image is repeated across 3 channels

## CNN Architecture:

Model: "sequential"

Layer (type)	Output Shape	Param #
model (Model)	(None, 7, 7, 512)	14714688
flatten (Flatten)	(None, 25088)	0
dropout (Dropout)	(None, 25088)	0
dense (Dense)	(None, 1024)	25691136
dropout_1 (Dropout)	(None, 1024)	0
dense_1 (Dense)	(None, 512)	524800
dropout_2 (Dropout)	(None, 512)	0
dense_2 (Dense)	(None, 256)	131328
dropout_3 (Dropout)	(None, 256)	0
dense_3 (Dense)	(None, 1)	257
Total params: 41,062,209		
Trainable params: 28,707,329		
Non-trainable params: 12,354,880		

Model info:





- The model contains out of 2 parts:
  - Part 1) A pre-trained [VGG16 model](#):
    - A 16 layer pre-trained convolutional neural network based on imagenet. The weights are not trainable.
  - Part 2) Four dense layers, each with its own dropout layer

## 3. Algorithm Training

### Parameters:

- Types of augmentation used during training
  - Rescale: pixel\_value / 255
  - Horizontal flips: True,
  - Vertical flip: False,
  - Height shift range: 0.1,
  - Width shift range: 0.1,
  - Rotation range: 20,
  - Shear range: 0.1,
  - Zoom range: 0.1

- Batch size: 64
- Optimizer: Adam
- Optimizer learning rate: 0.0003
- Layers of pre-existing architecture that were frozen: First 16 layers of a VGG model
- Layers of pre-existing architecture that were fine-tuned: None
- Layers added to pre-existing architecture: Flatten, Dense and Dropout layers



### Final Threshold and Explanation:

- Threshold: 0.37415
- F1 Score: 0.46824

The final threshold of 0.37415 was based on the highest F1-Score of 0.46824. Based on [this paper: CheXNet: Radiologist-Level Pneumonia Detection on Chest X-Rays with Deep Learning](#), the average radiologist had a F1 Score of 0.387. This means that our model achieved a better F-1 score than the average radiologist. *(Please note that we have no information about who these radiologists were and whether they are representative of all radiologists..)*

## 4. Databases

- Dataset: [NIH Chest X-ray Dataset](#)
- Amount of Chest, X-ray images: 112,120
- Meta data:
  - Image Index
  - Finding Labels
  - Follow-up
  - Patient ID
  - Patient Age
  - Patient Gender
  - View Position
  - Original Image Size
  - Original Image Pixel Spacing
- Amount of images with Pneumonia: 1431 (1.28%)
- Amount of images without Pneumonia: 110689 (98.72%)

### Description of Training Dataset:

1145 images (80%) which contained Pneumonia will be used in the training dataset.

- Amount of Chest, X-ray images: 2290
  - Images containing pneumonia: 1145 (50%)
  - Images without pneumonia: 1145 (50%)

**Description of Validation Dataset:** The remaining 286 images (20%) which contained Pneumonia will be used in the validation dataset.

- Images containing pneumonia: 286 (25%)
- Images without pneumonia: 858 (75%)

## 5. Ground Truth

- ChestX-ray dataset comprises 112,120 frontal-view X-ray images of 30,805 unique patients with the text-mined fourteen disease image labels (where each image can have multilabels)
- The disease labels were created using Natural Language Processing (NLP) to mine the associated radiological reports
- The labels are expected to be >90% accurate and suitable for weakly-supervised learning
- The data includes 14 common thoracic pathologies:
  - Atelectasis
  - Consolidation
  - Infiltration
  - Pneumothorax

- Edema
- Emphysema
- Fibrosis
- Effusion
- Pneumonia
- Pleural thickening
- Cardiomegaly
- Nodule
- Mass
- Hernia

For more info: <https://www.kaggle.com/nih-chest-xrays/data>

## 6. FDA Validation Plan

### Patient Population Description for FDA Validation Dataset:

- The patient can be men or women
- The recommended age is between 1 and 90 years old
- Chest X-Ray image must be taken in the AP or PA position
- Chest X-Ray image must be in DICOM format

### Ground Truth Acquisition Methodology:

- Silver Standard: Validation by 4 different radiologists

### Algorithm Performance Standard:

- The algorithm's F1 score should be more than that of average radiologist (0.387)